Affinity Selection Mass Spectrometry (ASMS) High-Throughput Screening of the GLP1 Receptor Zane Thistleford¹, Patrick O'Loughlin¹, Donald Green¹, Arrin Katz¹, Jan Kubicek², Barabara

Maertens² and Can Özbal¹

¹Momentum Biotechnologies, 3 Federal St., Suite 300, Billerica, MA 01821 ²Cube Biotech, Creative Campus Allee 12, 40789 Monheim, Germany

Abstract

High-throughput screening of membrane-bound proteins can be particularly challenging, because hydrophobic test compounds bind non-specifically to the artificial membranes and detergent micelles typically used to solubilize and stabilize target constructs. In this study, we describe a successful highthroughput screen against the GLP1 receptor (GLP1R), an important target in the treatment of type 2 diabetes and obesity. The NativeMP[™] platform was used to prepare the target construct, which consisted of the full-length GLP1R protein expressed in a copolymer-based nanodisc. A customized ASMS workflow was developed to enable compounds to be screened simultaneously against the GLP1R nanodisc (target screening) and a copolymer-matched blank nanodisc (off-target screening). Customized hit selection criteria were implemented using Aston ASMS data analysis software, and test compounds that displayed higher affinity for GLP1R than for the blank nanodisc were selected as hits. This screening configuration enabled test compounds that non-selectively bind the nanodisc construct to be efficiently eliminated. A successful pilot project (4K compounds) identified 5 unique stoichiometric binders, while a subsequent diversity screen (50K compounds) identified 75 hits (0.15% hit rate), of which 34 were confirmed in secondary assays. This screening approach can be applied to a wide range of membrane-bound protein targets, including GPCRs and ion channels.

Advantages of ASMS

Minimal assay development required	Ideal platform for selectivity screening					
High throughput	100K compounds screened in pools of 320 in 38 hours					
Low protein / RNA usage	 Primary screening (100K compounds in pools of 320) requires less than 1mg of a 50 kD target protein. 					
Solution-phase binding with native target	No need for target tagging or surface immobilization					
Ability to work with large targets (including protein/protein and protein/oligonucleotide complexes)	Large targets can be challenging for SPR or interferometry					
Hits are immediately available for confirmation and follow up	 Unlike DEL screening, the identity of hits are known, and hit compounds can be available as powders 					
Signal level is quantitatively correlated to target occupancy	Subtle changes in binding affinity can be detected in selectivity screening					
HTS Project Workflow Chemical Diversity Library 16 pools of 250 test compounds for pilot 200 pools of 250 test compounds for HTS Off-Target Screen 0.25 mM Target Protein Off-Target Screen 20 mM Blank Nanodisc						
Acquire ASMS data in parallel						
Identify test compounds unique to Target						
Dose-response experiments						
and Kd						
Screening and Data Analysis Parameters						
ASMS Screening Conditions Hit (Criteria Applied for Aston Software					
Sample volume 10 µL Peak Area	Min 5x background					

Cross-well filter

Off-target filter

Score

Blank Subtraction Min 4x peak area for target vs DMSO

Min 2x peak area target vs off-target

Max 3

Min 60

0.25 μ**Μ**

250

GLP1 Target Conc

Pooling Density

Blank Nanodisc Conc $20 \mu M$

Test Compound Conc 1 μ**M** each

GLP1R – Type 2 Diabetes Target



Some GLP1R therapies reduce the risk of heart attack and stroke, which is valuable for patients cardiovascular disease

Cardiovascular Protection

Oral Availability

Small molecules can be taken orally, eliminating the need for injections and making treatment for diseases like type 2 diabetes and obesity more convenient



NativeMPTM Platform Stabilization of GLP1R

- Full-length
- First ever stabilized wild-type
- Retention of native host cell lipid environment
- Function and stability preserved
- Available in mg amounts
- Compatible with assay development at body temperature

Sample Purity SDS-PAGE kDa 95 72 ____ _____



Unique characteristics of Lipid-Membrane Protein Particles

Detergentstabilized Membrane Protein

- Non-lipid environment • Function and stability
- compromised **Restricted** assay
- temperatures
- High-risk project with long timeline
- Not compatible with full-length wild-type GPCRs



NativeMP[™] platform

Lipid Membrane Protein Particle in Synthetic Polymers

- native host cell lipid environment
- **Function and stability** preserved
- Assays can be performed at body temperature
- Lower risk and faster project timeline
- **Enables analysis of** full-length wild-type GPCRs

Homogenous sample without aggregation



DLS-analysis





Hit Confirmation Assay

Example screenshot of custom Aston ASMS Software for data analysis. Each test compound is

MS signal quality, and signal must be unique to the the target well. Graphical displays of each

interrogated in the target well in parallel with DMSO-only controls, matched off-target injections (in

this case, blank nanodisc), and non-target wells. Primary hits must pass user-selected criteria for

compound's Extracted Ion Chromatogram and Mass Spectrum, as well as each injection's SEC-UV

0.00 0.00 0.00 0.00 0.00 0.00 0.00 0.00 0.00 0.00 0.00 0.00

0.00 0.00 00.0 00.0

0.00 0.00 0.00 0.00 0.00 0.00 0.00 0.00

0.00 0.00 0.00 0.00

0.00 0.00 0.00

0.00 0.00 0.00 0.00



trace, are displayed to enable rapid manual QC of results as needed

Example results from a singleton hit confirmation assay using Agilent Mass Hunter Software. The **Extracted Ion Chromatogram** (Panel A) and Mass Spectrum (Panel B) of this test compound indicate that the structure has been detected and quantified. ASMS signal intensity is defined as the area under the curve for the compound chromatogram peak (shaded in green in Panel A) and corresponds to the concentration of the test compound. Test compound concentration is proportional to target occupancy, which is in turn proportional to the affinity of the ligand for the target.

1 2 3

Time (minutes)

Blank Method: PPM: 10 Abundance Threshold: 250 Nominal Zero: 0.1 Min Area: 2500 Min T: 1 Max T: 4

4

0.00 0.00 0.00

0.00 0.00 0.00

③ offtarget-E12 0.00

Na: Peak area too low 113.41 Peak area ratio too lo

Results

Initial testing was performed using a set of 4,000 test compounds (16 pools x 250 compounds) selected at random from the Momentum diversity library (300,000 compounds). Assays were run using a range of GLP1R and blank nanodisc concentrations to optimize parameters. Ultimately, a target concentration of 0.25 uM and an off-target (blank nanodisc) concentration of 20 uM were selected, resulting in 6 verified hits (0.15% hit rate). These conditions were then used to perform a primary ASMS screen with 49,546 test compounds (Momentum PS7 sublibrary). This library consists of 199 mass-coded pools of 250 test compounds each, assorted such that each test compound in a given pool has a mass difference ≥ 0.2 g/mol from all pool-mates. The primary screen was completed in under 24 hours and identified 75 hits (0.15% hit rate), which were further interrogated in a confirmation assay. Primary hits were run as singletons against the target and blank nanodisc, resulting in 34 confirmed hits for an overall confirmation rate of 45%. Confirmed hits were subjected to 8-point dose-response experiments using serial dilutions starting at 50 uM and including DMSO-only and no-target controls. Many test compounds exhibited saturating dose-response curves, indicating specific and stoichiometric binding to the target protein.

Target Selectivity of Hits

Cnd ID	Molecular Formula	Target ASIVIS Peak	Off-target ASIVIS	Ratio: Target/Off-
еритр	inforceular ronnata	Area	Peak Area	target
1271517	C19 H17 N3 O	119486	ND	NA
1271520	C18 H14 Cl N3 O	81644	ND	NA
1272547	C19 H23 N3 O	37321	ND	NA
1275201	C19 H31 N5	555608	ND	NA
1275764	C14 H12 N2 O S2	18017	ND	NA
1275958	C16 H21 N3 O3 S	7797	ND	NA
1290789	C20 H19 N3 O3	87402	ND	NA
1291652	C21 H30 F N3 O	133181	ND	NA
1292203	C21 H27 N5	85755	ND	NA
1293841	C14 H13 Cl N4 O	32335	ND	NA
1299453	C24 H29 N3 O2	96432	ND	NA
1301285	C23 H24 N4 O4	21947	ND	NA
1305682	C21 H25 N5 O3	116767	ND	NA
1306039	C26 H28 F N3 O2	401944	ND	NA
1306554	C22 H27 N5 O S	44319	ND	NA
1310492	C18 H30 N2 O3 S	47378	ND	NA
1310653	C19 H25 N3 O2	128147	ND	NA
1311222	C19 H21 Cl N2 O2	24869	ND	NA
1311531	C16 H20 F N3 O S	44724	ND	NA
1311680	C16 H22 N4 O S	228963	ND	NA
1312044	C22 H23 N3 O2	276813	ND	NA
1314345	C18 H21 N7 S	224715	ND	NA
1314986	C19 H17 F N4 O	16347	ND	NA
1315345	C22 H27 N5	49095	ND	NA
1306217	C27 H27 F N4 O	639956	28490	22.5
1291709	C17 H28 N6 O	621132	63855	9.7
1274598	C19 H24 F3 N5 O S	76527	10939	7.0
1293685	C22 H24 F N O3	80234	14928	5.4
1306507	C22 H32 N6 O2	894676	167319	5.3
1272721	C16 H15 F N4 O S	340691	93393	3.6
1290335	C17 H18 Cl N5	65722	21021	3.1
1309913	C14 H16 F2 N2 O3	26181	8574	3.1
1277620	C18 H20 N6 O S	51366	17329	3.0
1298002	C22 H23 N5 O	72438	27082	2.7

Hit confirmation assays were run against both 0.25 uM GLP1R target protein and 20 uM blank nanodisc. The peak area for each injection is shown in the table to the left. Of 34 confirmed hits, 24 bound only the target protein, with no detectable binding to the blank nanodisc. The remaining 10 confirmed hits bound both the target protein and blank nanodisc to some extent, suggesting weak non-specific binding to the nanodisc lipids in addition to specific target interactions. A target peak area of at least 2X the off-target peak area was required for hit selection.



Affinity Determination

Eight-point serial dose-response experiments were run on the 34 confirmed hits to determine the binding affinity of each test compound. A saturating binding curve indicates specific and stoichiometric binding to the target, from which the binding affinity can be determined. A linear dose-response curve (as seen in Panel C) indicates super-stoichiometric binding, likely due to nonspecific interactions between the target and test compound; this compound represents a true binder but is unlikely to be a useful lead compound. Kd values for the 5 specific binders ranged from 1.7 uM to 45 uM. Interestingly, two test compounds (shown in Panels B and F) share over 97% structural homology but exhibit a three-fold difference in Kd values, providing early insight into the structure-activity relationship.

Conclusions

Membrane-bound proteins pose unique challenges for binding assays due to extensive nonspecific binding between poorly soluble test compounds and lipids required to stabilize the target protein. We developed a rapid, robust and automated binding assay that overcomes these challenges to target a novel GLP1R-nanodisc construct. The use of the nanodisc system minimizes the lipid load used in the assay, while the off-target screening approach effectively eliminates nonspecific binders from the hit pool, greatly streamlining downstream interrogation and increasing the percentage of compounds useful for a drug discovery campaign.

www.momentum.bio www.cube-biotech.com