



Sulfo-DIBMA



Product	Catalog No.	Package size
Sulfo-DIBMA (10 x 50 mg)	18531	10 x 50 mg
Sulfo-DIBMA (1 g)	18532	1 g
Sulfo-DIBMA (10 x 1 g)	18533	10 x 1 g
Sulfo-DIBMA (50 g)	18534	50 g



Product Description

The use of a Poly(diisobutylene-alt-N',N'-dimethyl(maleimidopropyl))ammonium propane sulfonate copolymer (Sulfo-DIBMA) for stabilization of membrane proteins can provide bicelles with membrane proteins from native membranes in the absence of detergents. This is achieved by wrapping around a patch of a lipid bilayer to form a disc-like particle or nanodisc. The Sulfo-DIBMA-based products contain the copolymer and a 50 mM HEPES buffer, adjusted to pH 7.5, so only dd water has to be added for direct application. The pH value has been selected to be very effective for protein solubilization.

Cube Biotech's Sulfo-DIBMA is a highly purified electroneutral copolymer, with a molecular weight (Mw) of ~8.000. After dissolving the lyophilized copolymer powder with membrane protein-containing buffer, the concentration should range from 1.0 to 5.0%. Copolymers provide a hydrophobic surface facing the lipids and a hydrophilic surface on the outside. This setup makes nanodiscs highly soluble in aqueous solutions and allows the solubilization of membrane proteins in the absence of detergents. This product can be used with phospholipids, such as dimyristoyl-glycero-phosphocholine (DMPC) or palmitoyl-oleoyl-phosphatidyl-choline (POPC) in combination with sodium cholate.

The complex from Sulfo-DIBMA and membrane protein can be used in many biophysical assays, such as SDS-PAGE, SEC, Western Blot, UV/Vis spectroscopy, and many chromatographic procedures.

Reconstitution of the copolymer solution:

Sulfo-DIBMA copolymers are delivered as a lyophilized solution containing 50 mM HEPES, pH 7.5. Each aliquot contains 50 mg of polymer, 1 g, 10 g or 50 g respectively. Adding 0.5 ml double distilled water per 50 mg of polymer will restore the original solution with a copolymer concentration of 10%. This stock can be diluted further as required by different application protocols.

Technical Details

Name	Poly(diisobutylene-alt- N^+ , N^- -dimethyl(maleimidopropyl))ammonium propane sulfonate copolymer, sodium salt in 50 mM HEPES, pH 7.5
Solubility	>10% (H ₂ O)
Color	White to slightly yellow
Odor	Odorless
pH (dissolved)	7.5 ± 0.3

Shipping & Storage

Shipping temperature	Ambient temperature
Long-term Storage (lyophilized copolymer)	-20°C for several years
Short-term Storage (dissolved copolymer)	2-8°C for several days

Additional Information

For Sulfo-DIBMA protocols, please visit our webpage at: <https://cube-biotech.com/products/membrane-protein-stabilization/synthetic-nanodisc-products/>

Disclaimer

Patent Pending

The product described above is based on a substance published within the work of D. Glueck et al., Small, 2022, 2202492, DOI: 10.1002/smll.202202492.

Our products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.