


HisCube Ni-INDIGO His-tag MagBeads MINI Kit

Kit contents

1 x 5 ml PureCube Ni-INDIGO MagBeads
2 x 50 ml Binding Buffer
2 x 50 ml Wash Buffer
1 x 50 ml Elution Buffer
1 x Cube MagBead Separator
1 x 500 mM DTT (add 1,016 ml ddH₂O)
1 x 500 mM EDTA (add 0.866 ml ddH₂O)

 Storage Temperature 4 °C / 39 °F

 Protocol



Scan to download
the full protocol

Fast Protocol

This „Fast Protocol“ is only meant as reminder for experienced HisCube Ni-INDIGO MINI Kit users

- 1 Pipet 200 μ l Ni-INDIGO resin slurry into your reaction tube of 1.5 or 2 ml volume.
- 2 Equilibrate the MagBeads with 1.5 ml Binding Buffer by shaking the reaction tube for about **30 seconds** end-over-end. **Repeat** this step. Separate with the magnet and discard the supernatant.
- 3 Immediately before loading: Re-filter the cleared lysate through a 0.2 μ m filter e.g. a syringe filter.
- 4 Load the required volume (max. 1.8 ml) of filtered cleared lysate. **Vortex for 15 sec.**
- 5 Repeat the vortexing of Step 7 every **15 minutes** for **1 hour** at **4 °C** or for **30 minutes** at room temperature.
- 6 Separate with the magnet for at least **20 sec.** Save the flow-through in a separate tube and place the reaction tube containing the purification MagBeads back into the magnet separator.
- 7 Load the tube with up to 1.2 ml of Wash Buffer and separate for at least **20 sec.** Remove the supernatant. **Repeat** at least two times. Save all supernatants.
- 8 Elute the target protein by adding up to 100 μ l Elution Buffer for **10 to 30 min** and then separate for at least **20 sec.** If necessary, repeat the elution step up to five times.