

PureCube Epoxy Activated Agarose



Product	Catalog No.	Package size
PureCube Epoxy Activated Agarose (10 mL)	50703	20 mL 50% suspension
PureCube Epoxy Activated Agarose (50 mL)	50705	100 mL 50% suspension
PureCube Epoxy Activated Agarose (250 mL)	50710	500 mL 50% suspension
PureCube Epoxy Activated Agarose (500 mL)	50712	1000 mL 50% suspension

Product Description

PureCube Epoxy Activated Agarose has been synthesized for the direct covalent coating for affinity purification.

This affinity chromatography matrix is based on BioWorks Workbeads, consisting of 6% cross-linked agarose. The material is highly porous to allow for optimal protein interaction. Cross-linked agarose is also physically very stable, making it suitable for purification processes under low pressure with flow rates up to 6 mL/min (optimal 0.5 – 2 mL/min). Our agarose is very homogeneous in size with a medium particle diameter of 100 µm, yielding a high degree of reproducibility between individual purification runs.

An Epoxy Activated-modified epoxide function is coupled to the agarose beads with a C4 spacer to obtain a matrix with highest binding capacity for amino and thiol functions. The epoxy group density is higher than 20 µmol/mL, as determined by acidimetric titration. PureCube Epoxy Activated Agarose is delivered as a 50% suspension. Therefore, 1 mL suspension will yield a 500 µL bed volume. The suspension contains 100% isopropanol to prevent hydrolysis and microbial growth.

Shipping & Storage

Shipping Temperature	Ambient temperature
Short-term Storage	In equilibration buffer (see protocol) at 4 °C
Long-term Storage	Best lyophilized, or in 100% isopropanol at 4 °C