

## Sulfo-Cubipol Lite



Product	Catalog No.	Package size
Sulfo-Cubipol Lite (10x50mg)	19001	5 g
Sulfo-Cubipol Lite (1g)	19002	1 g
Sulfo-Cubipol Lite (10x1g)	19003	10 g
Sulfo-Cubipol Lite (50g)	19004	50 g

### Product Description

The use of a Sulfo-Cubipol Lite HEPES for stabilization of membrane proteins could provide bicelles with membrane proteins from native membranes in absence of detergents, by wrapping around a patch of a lipid bilayer to form a disc-like particle or nanodisc. The Sulfo-Cubipol Lite HEPES based products contain the copolymer and a 50 mM HEPES buffer, adjusted to pH 7.5, so only dd water has to be added for direct application. The pH value has been selected being very effective for protein solubilization. Sulfo-Cubipol Lite HEPES from Cube Biotech is a highly purified copolymer, with a molecular weight (Mw) of ~ 9000 Da. After dissolving the lyophilized copolymer powder with membrane protein-containing buffer, the concentration should be in the range from 1.0 to 5.0%. Copolymers provide a hydrophobic surface facing the lipids, and a hydrophilic surface at the outside. This setup makes nanodiscs highly soluble in aqueous solutions and allows the solubilization of membrane proteins in the absence of detergents. The product can be used with phospholipids, such as dimyristoyl-glycero-phosphocholine (DMPC) or palmitoyl-oleoyl-phosphatidyl-choline (POPC) in combination with sodium cholate.

The complex from Sulfo-Cubipol Lite HEPES and membrane protein can be used with many biophysical assays, such as SDS-PAGE, SEC, Western Blot, UV/Vis spectroscopy, and many chromatographic procedures.

### Reconstitution of Copolymer Solution

Sulfo-Cubipol Lite HEPES copolymers, as delivered, are lyophilized from a solution containing 50 mM HEPES, pH 7.5. Each aliquot contains 50 mg of copolymer, 1 g, 10 g, or 50 g, respectively. Adding 0.5 ml double-distilled water per 50 mg of copolymer will restore the original solution with a copolymer concentration of 10%. This stock can be diluted further as required by the different application protocols.

**Technical Details**

Name	Sulfo-Cubipol Lite, sodium salt in 50 mM HEPES, pH 7.5
Solubility	>10% (H <sub>2</sub> O)
Color	White to slightly yellow
pH (dissolved)	7.5 ± 0.3

**Shipping & Storage**

Shipping Temperature	Ambient Temperature
Storage of lyophilized copolymer	-20°C for several years
Storage of dissolved copolymer	2-8°C for several days

**Additional Information**

For Sulfo-Cubipol Lite protocols, please visit our corresponding product websites under: <https://cube-biotech.com/products/membrane-protein-stabilization/copolymer-nanodisc-products/cubipol/>

**Disclaimer**

Our products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.