

Coupling of proteins to PureCube Carboxy Agarose

Overview

This protocol delineates a coupling procedure for proteins to PureCube Carboxy Agarose using the cross-linking reagents EDC and NHS. Proteins are coupled covalently and can be used for different applications, e.g. purification of interaction partners.

Amounts given in this protocol are for 2 mL of a 50% Carboxy Agarose suspension, corresponding to 1 mL agarose bed volume. The process can be linearly scaled up and down from 100 μ L to several 100 mL bed volume.

Separation of Agarose and supernatants in low mL scale can be done by centrifugation, or, more conveniently, using PureCube 1-step batch Columns. See the separate protocol available for a combination of these products. For larger scales, use of special equipment like glass vacuum tank filters, is recommended.

Please contact us if you have questions or need assistance optimizing a protocol for your application (contact@www.cube-biotech.com); other protocols can also be found at www.cube-biotech.com/protocols.

Equipment

- Centrifuge for 15 mL tubes
- Centrifuge tubes (15 mL)
- End-over-end mixer or thermomixer
- Spectrophotometer

Materials

- PureCube Carboxy Agarose (10 mL, Cube Biotech #50103)
- Sodium dihydrogen phosphate
- Sodium chloride
- Sodium acetate trihydrate
- Sodium hydroxide (NaOH)
- Acetic acid
- Ethanol
- Ethanolamine
- N-(3-Dimethylaminopropyl)-N'-ethyl-carbodiimide hydrochloride (EDC)
- N-Hydroxysuccinimide (NHS)

Solutions and buffers

Phosphate buffer, pH 6.0, 250 mL

| Component | Final concentration | Molecular weight (g/mol) | Stock concentration | Amount needed for buffer |
|--|---------------------|--------------------------|---------------------|--------------------------|
| NaH ₂ PO ₄ dihydrate | 25 mM | 156.01 | n.a. | 975 mg |
| Instructions: Dissolve sodium phosphate in 200 mL water, adjust the pH to 6.0 with NaOH. Add water to a total volume of 250 mL. | | | | |

PBS Buffer, pH 7.2, 250 mL

| Component | Final concentration | Molecular weight (g/mol) | Stock concentration | Amount needed for buffer |
|--|---------------------|--------------------------|---------------------|--------------------------|
| NaH ₂ PO ₄ dihydrate | 150 mM | 156.01 | n.a. | 5.85 g |
| NaCl | 100 mM | 58.44 | n.a. | 1.463 g |
| Instructions: Dissolve components in 200 mL water, adjust the pH to 7.2 with NaOH. Add water to a total volume of 250 mL. | | | | |

Quenching Buffer, pH 7.4, 250 mL

| Component | Final concentration | Molecular weight (g/mol) | Stock concentration | Amount needed for buffer |
|--|---------------------|--------------------------|---------------------|--------------------------|
| Ethanolamine | 1 M | 61.08 | n.a. | 15.27 g |
| Instructions: Dissolve component in 200 mL water, adjust the pH to 7.4 with HCl. Add water to a total volume of 250 mL. | | | | |

Agarose Storage Buffer, pH 6.5, 250 mL

| Component | Final concentration | Molecular weight (g/mol) | Stock concentration | Amount needed for buffer |
|---|---------------------|--------------------------|---------------------|--------------------------|
| Sodium acetate trihydrate | 20 mM | 136.08 | n.a. | 674 mg |
| Ethanol | 20 % (v/v) | | 100 % (v/v) | 51 mL |
| Instructions: Dissolve sodium acetate in 150 mL water, adjust the pH to 6.5 with acetic acid. Add 48 mL water and 51 mL ethanol to yield a total volume of 250 mL. | | | | |

Procedure

1. Transfer 2 mL PureCube Carboxy Agarose suspension (corresponding to 1 mL bed volume) into a 15 mL centrifuge tube.
2. Spin the tube at 500 x g to pellet the agarose. Remove the supernatant. Resuspend the agarose in 3 mL Phosphate Buffer.
3. Dissolve 250 mg NHS in 1 mL Phosphate Buffer, add it to the agarose suspension, and mix by vortexing.
4. Dissolve 250 mg EDC in 1 mL Phosphate Buffer, add it to the agarose suspension, and mix by vortexing.
5. Incubate at room temperature for 1 h on an end-over-end shaker or in a thermoshaker.
6. Prepare a solution of 2.5 mL PBS containing the protein to be coupled to the agarose. The exact protein amount needs to be optimized, and 5 to 15 mg protein is a good starting point.
7. Add the protein solution to the agarose and mix by vortexing.
8. Depending on the temperature stability of the protein, incubate at room temperature or 4°C for 2 h on an end-over-end shaker or thermoshaker.
9. Spin the tube at 500 x g to pellet the agarose. Remove the supernatant and analyze the supernatant in a spectrophotometer. Record absorbance at 280 nm to monitor coupling efficiency.
10. Add 5 mL double distilled water to the agarose pellet, mix by vortexing, and spin at 500 x g. Remove the supernatant.
11. Repeat step 10 five times.
12. Add 5 mL Quenching Buffer and incubate again for 1 h at room temperature or for 4 hours at 4°C.
13. Wash four times with 5 mL PBS each, and twice with 5 mL double distilled water each.
14. Resuspend the coupled Agarose in 2 mL Agarose Storage buffer, yielding a 50% suspension. Store at 4°C.

Tip: The coupling reaction can be linearly scaled up and down, by increasing or decreasing the amounts of buffers and solutions described in this protocol.

Note: NHS and EDC should always be prepared fresh. Equilibrate the two chemicals to room temperature before weighing, and store the powders under protective gas (nitrogen). Add the two chemicals immediately one after another to prevent hydrolysis of the NHS-activated matrix.

Tip: When coupling a particular protein for the first time, try 3-5 different protein concentrations to make sure you are offering enough protein in the reaction but not wasting any protein.

Tip: Monitoring absorbance at 280 nm tells you about the coupling efficiency of the protein (compare A₂₈₀ of the original protein solution to the supernatant in step 6 to determine % coupling). It also helps you identify the optimal amount of protein required for efficient coupling.

Note: The quenching step ensures that no free NHS groups are left on the agarose matrix that might interfere with subsequent assays.

