

MSP1E3D1-His “Ultra Bright”, human

Cat. No.: 26652

Description	MSP1E3D1-His “Ultra Bright” is a His-tagged membrane scaffold protein (MSP) derived from apolipoprotein A-I that wraps around a phospholipid bilayer patch to form stable, water-soluble nanodiscs of 12–14 nm in diameter, enabling the stabilization of membrane proteins in a defined lipid environment. Its defining feature is the addition of three tryptophan residues in sequence, resulting in a significantly increased intrinsic fluorescence of the scaffold protein itself. This product contains a His-tag, keep that in mind when choosing the affinity tag for your membrane protein.
Form	Lyophilized
Amount	2 mg
Available origins	Human
Purity	> 90% (Determined by SDS-PAGE)
Number of amino acids	277
Molecular mass	32.81 kDa
Extinction coefficient (in H₂O) ϵ_{280}	Ext. coefficient 50420 Abs 0.1% (=1 g/l) 1.537
Buffer	20 mM TRIS pH 7.4, 100 mM NaCl, 0.5 mM EDTA
Affinity tag	His tag
Sequence	GHHHHHHH DYDIPTT ENLYFQG STFSKLREQLG PVTQEFWDNLEKETEGLRQEMS KDLEEVKAKVQ PYLDDFQKKWQEEMELYRQKVE PLRAELQEGARQKLHEL WE KLS PLGEEMRDRARAHVDALRTHLA PYLDDFQKKWQEEMELYRQKVE PLRAELQEGARQKLHEL WE KLS PLGEEMRDRARAHVDALRTHLA PYSDELQRQLAARLEA W KENG ARLA EYHAKATEHLSTLSEKAK PALEDLRQGLL PVLESFKVSFLSALEE W TKKLNTQ
Stability & Storage	Lyophilized protein: stable at -20 °C for 24 months. Reconstituted protein: stable at 2-8 °C for several days.
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP].

	A Material Safety Data Sheet is provided.
Application	The increased intrinsic fluorescence of MSP1E3D1-His "Ultra Bright" allows easy detection and monitoring of nanodisc properties without exogenous labels, and enables its use as an intrinsic FRET donor for labeled proteins. It further supports binding studies via analytical ultracentrifugation (AUC) and microscale thermophoresis (MST) without the need for extrinsic fluorophores. Beyond its fluorescence advantages, MSP1E3D1-His "Ultra Bright" is suited for Cryo-EM, ELISA, protein-protein interaction studies, and cell-free expression systems, supporting a broad range of biomedical, proteomic, and protein biochemistry applications.
Recommended usage	ELISA: dilute 1:100
Manufacturer	Cube Biotech

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.