

PureCube HiCap Strep-Tactin® FPLC column

Cat. No.: 34302, 34304, 34306, 34308

Description	PureCube HiCap Strep-Tactin® FPLC columns are prepacked with PureCube HiCap Strep-Tactin® agarose. The FPLC columns are applicable for purification of Strep-tag®II or Twin-Strep-tag® fusion proteins via chromatography workstations using 10-32 connections, syringes, or peristaltic pumps.			
Support	4% agarose, highly cross-linked			
Form	Pre-packed in PBS (100 mM sodium phosphate, 150 mM NaCl, pH 7.2) with 0.02% sodium azide			
Bead size	50-150 µm, spherical			
Binding capacity	5 mg protein/ml resin. Binding capacity was determined with a 82.6 kDa CD56 antigen. Please note that the binding capacity is protein dependent.			
Compatible reagents	2% Tween 20, 0.5% IGEPAL® Ca-630/Nonidet P40, 1% Dodecyl-β-D-maltoside, 1 mM PMSF, 10% ethanol, 5 M NaCl, 2 M (NH ₄) ₂ SO ₄ , 22% glycerol.			
Application	The PureCube FPLC column is compatible with almost all commonly used FPLC systems like ÄKTA™.			
	FPLC column volume	Physical dimension	Flow rate	Pressure stability
	1 ml (34302, 34304)	Inner diameter: 0.7 cm Height: 2.5 cm	1-3 ml/min (max. 6 ml/ml)	Maximum 3.5 bar
	5 ml (34306, 34308)	Inner diameter: 1.6 cm Height: 2.5 cm	1-2 ml/min (max. 6 ml/ml)	Maximum 3.5 bar
Stability	6 months after shipping			
Storage	2-8 °C			
Shipping	Room temperature			
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.			
Manufacturer	Cube Biotech			

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at contact@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.