

Strep Antibody HRP

Product	Catalog No.	Package size
Strep Antibody HRP (75 µl)	40071	1 x 75 µl

Product Description

75µl Anti-Strep-tag® II mouse monoclonal antibody, purified and conjugated to horseradish peroxidase.

Lot #	40071-15090082
Antibody Name	Strep Antibody HRP
Subclass	IgG
Form	Suspension in PBS containing 5 mg/ml bovine serum albumin; approx. 0.05% ProClin 300 and stabilizing agents
Specificity	Strep-tag® II; the antibody was produced in mouse using Strep-tag® II peptide conjugated to keyhole limpet hemocyanin as immunogen
Detections	N- or C-terminal Strep-tag® II and Twin-Strep-tag®
Dilutions	Dilute 1:10000 or higher for detection of Strep-tag® II fusion proteins coated to microplates. Dilute 1:4000 for chromogenic and 1:32000 for chemiluminescence detection of Strep-tag® fusion proteins after Western blotting
Amount	0.075 ml

Shipping & Storage

Shipment Temperature	Shipped at room temperature (RT)
Storage	Recommended: 2-8°C

Disclaimer: Our products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

Strep Antibody is manufactured by IBA GmbH.

Trademarks: Strep-tag® (IBA GmbH)