

PureCube HiCap Strep-Tactin® MagBeads XL

Cat. No.: 57405, 57425

Description	PureCube HiCap Strep-Tactin® MagBeads XL are magnetic beads with a diameter of 90 µm (average) coupled with Strep-Tactin®. They can be used for affinity purification of proteins with the Strep-tag®II or Twin-Strep-tag®.
Support	Spherical magnetic agarose beads, cross-linked, 6% agarose
Form	5% (v/v) suspension in PBS (100 mM sodium phosphate, 150 mM NaCl, pH 7.2) with 0.02% sodium azide (1000 µL suspension corresponds to 50 µL magnetic beads)
Bead size	71 – 112 µm (90 µm average)
Binding capacity	7 mg protein/ml beads. Binding capacity was determined with a 48 kDa Kinase. Please note that the binding capacity is protein dependent.
Compatible reagents	2% Tween 20, 0.5% IGEPAL® Ca-630/Nonidet P40, 1% Dodecyl-β-D-maltoside, 1 mM PMSF, 10% ethanol, 5 M NaCl, 2 M (NH ₄) ₂ SO ₄ , 22% glycerol.
Application	Elution can be performed by boiling in denaturing SDS gel loading buffer or by application of desthiobiotin containing buffer.
Stability	12 months after shipping
Storage	2-8 °C
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Manufacturer	Cube Biotech

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.