

MSP Starter Kit for In Vitro Protein Synthesis (-His)

Cat. No.: 26702 (8x7.5 µL)

Description	The MSP Starter Kit for In Vitro Protein Synthesis is ready to use for cell-free protein synthesis: it contains 8 different pre-assembled MSP-nanodiscs, which can be added directly to the expression reaction so the membrane protein inserts co-translationally into a membrane with predefined lipid composition, no detergents or post-synthesis reconstitution needed. This selection of MSP pre-assembled nanodiscs is particularly well suited for GPCRs. Two scaffold sizes (MSP1D1, MSP1E3D1) × four lipid compositions (DOPG; DOPG/Cardiolipin; POPC/POPS; POPC/POPS/Cholesterol) allow quick identification of the best-performing condition, manually or on a liquid handler. The His-tag enables fast affinity purification, and matching Scale-up vials carry the winning condition forward to µg-to-mg quantities of purified, functional membrane protein.
Form	pre-assembled MSP nanodisc
Amount	8 x 7.5 µL (500 µM)
Available origins	Human
Purity	> 90% (Determined by SDS-PAGE)
Number of amino acids	MSP1D1 - His: 217 MSPE3D1-His: 277
Molecular mass	MSP1D1-His: 25.3 kDa MSPE3D1-His: 32.6 kDa
Extinction coefficient (in H₂O) ε₂₈₀	MSP1D1-His: 25,440 M ⁻¹ cm ⁻¹ MSPE3D1-His: 29,400 M ⁻¹ cm ⁻¹
Buffer	20 mM TRIS pH 7.4; 100 mM NaCl; 0.5 mM EDTA
Affinity tag	N-Terminal His-tag
Stability & Storage	-80° C for several months Avoid repeated freeze-thaw cycles
Shipping	Dry ice
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

Application	Primarily used in cell-free (in vitro) synthesis to screen and identify the best-performing nanodisc condition for a given GPCR. The His-tag then enables fast, affinity-based purification of the reconstituted receptor, supporting downstream structural studies (e.g. cryo-EM), ligand-binding and functional assays, and early-stage drug discovery screening.
Recommended usage	In Vitro Protein Synthesis 5–75 μ M (0.5–7.5 μ L)
Manufacturer	Cube Biotech

For research use only**Trademark information**

The owners of trademarks marked by “®” or “TM” are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.