

GentleLys – Dissolve

Native Cell Lysis Buffer

Cat. No.: 18803, 18806

Description	GentleLys Dissolve is a 100% detergent-free, non-denaturing solution to lyse cultured insect and mammalian cells while solubilizing all membranes in the process.
Form	Lyophilized (1 ml) or solution (100 ml)
Application	The buffer works for eukaryotic cells (except for intact plant cells), not for bacterial cells. It is compatible with applications such as immunoassays, protein assays, mass spectrometry, and protein purification.
Stability	6 months after shipping
Storage	2-8 °C
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Manufacturer	Cube Biotech

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.