

MSP1D1-His “Ultra Bright”, human

Cat. No.: 26612 (2 mg)

Description	MSP1D1-His "Ultra Bright" is a his-tagged membrane scaffold protein (MSP) derived from apolipoprotein A-I that wraps around a phospholipid bilayer patch to form stable, water-soluble nanodiscs, enabling detergent-free solubilization and stabilization of membrane proteins in a near-native lipid environment. Its defining feature is the addition of three tryptophan residues in sequence, resulting in a significantly increased intrinsic fluorescence of the scaffold protein itself.
Form	Lyophilized
Amount	2 mg
Available origins	Human
Purity	> 90% (Determined by SDS-PAGE)
Number of amino acids	222
Molecular mass	26.1 kDa
Extinction coefficient (in H₂O) ϵ_{280}	Ext. coefficient 41940 Abs 0.1% (=1 g/l) 1.607
Buffer	20 mM TRIS pH 7.4, 100 mM NaCl, 0.5 mM EDTA
Affinity tag	His
Sequence	GHHHHHHH DYDIPTT ENLYF QG LKLLDNWDSVTSTFSKLREQL G PVTQEFWDNLEKETEGLRQEMS KDLEEVKAKVQ PYLDDFQKKWQEEMELYRQKVE PLRAELQEGARQKLHEL WE KLS PLGEEMRDRARAHVDALRTHLA PYSDELRLQRLAARLEA W KENG G ARLAEYHAKATEHLSTLSEKAK PALEDLRQGLL PVLESFKVSFLSALEE W TKKLNTQ
Stability & Storage	Lyophilized protein: stable at -20°C for several months. Reconstituted protein: stable at 2-8°C for several days.
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

Application	The increased intrinsic fluorescence of MSP1D1-His "Ultra Bright" allows easy detection and monitoring of nanodisc properties without exogenous labels, and enables its use as an intrinsic FRET donor for labeled proteins. It further supports binding studies via analytical ultracentrifugation (AUC) and microscale thermophoresis (MST) without the need for extrinsic fluorophores. Beyond its fluorescence advantages, MSP1D1-His "Ultra Bright" is suited for Cryo-EM, ELISA, protein-protein interaction studies, and cell-free expression systems, supporting a broad range of biomedical, proteomic, and protein biochemistry applications.
Recommended usage	ELISA: dilute 1:100
Manufacturer	Cube Biotech

For research use only**Trademark information**

The owners of trademarks marked by “®” or “TM” are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.