


REVISIONS			
ECOM	REV	DESCRIPTION	EFF DATE
1634	New	SL0107	9/24/10



Mat: PF 511 / PD 702

Tolerances unless otherwise specified:
 Without Decimal: ± 0.5 mm
 1 Place Decimal: ± 0.1 mm
 2 Place Decimal: ± 0.05 mm
 Material Thickness: ± 0.02 mm
 Undefined Radius: R0.2 mm
 Undefined Draft: 2°
 Vol: 2.2ml
 .3 depth excluded

All dims given from base of draft

PROPRIETARY AND CONFIDENTIAL		UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MM TOLERANCES ARE: ANGLES: $\pm 1^\circ$ SURF: ± 0.05 FRAC: ± 0.05	
THIS DRAWING IS THE PROPERTY OF AXYGEN SCIENTIFIC CORPORATION, IRVINE, CA, AND MAY NOT BE REPRODUCED, COPIED, OR DISCLOSED TO OTHERS WITHOUT WRITTEN CONSENT OF THE COMPANY. IT IS SUBJECT TO RETURN TO AXYGEN AND MAY NOT BE USED IN ANY WAY WITHOUT THE WRITTEN PERMISSION OF THE COMPANY.		APPROVALS	DATE
		DRAWN CRN	01/07/02
		CHECKED DT	01/07/02
		ENGR	
e:\drga\app\adrg\sl0107AA		DWG TITLE	P-2ML-SQ
SCALE 1:1, 5:1		SIZE	D
		DWG NO.	SL0107
		SHEET	1 OF 1