

MSP1E3D1 “Ultra Bright”, human

Cat. No.: 26662 (2 mg)

Description	MSP1E3D1 "Ultra Bright" is a membrane scaffold protein (MSP) derived from apolipoprotein A-I that wraps around a phospholipid bilayer patch to form stable, water-soluble nanodiscs, enabling detergent-free solubilization and stabilization of membrane proteins in a near-native lipid environment. Its defining feature is the addition of three tryptophan residues in sequence, resulting in a significantly increased intrinsic fluorescence of the scaffold protein itself.
Form	Lyophilized
Amount	2 mg
Available origins	Human
Purity	> 90% (Determined by SDS-PAGE)
Number of amino acids	256
Molecular mass	30.19 kDa
Extinction coefficient (in H₂O) ϵ_{280}	Ext. coefficient 47440 Abs 0.1% (=1 g/l) 1.571
Buffer	20 mM TRIS pH 7.4, 100 mM NaCl, 0.5 mM EDTA
Affinity tag	No tag
Sequence	G STFSKLREQLG PVTQEFFDNLEKETEGLRQEMS KDLEEVKAKVQ PFLDDFQKKFQEEMELFRQKVE PLRAELQEGARQKLHELQEKLS PLGEEMRDRARAHVDALRTHLA PFLDDFQKKFQEEMELFRQKVE PLRAELQEGARQKLHELQEKLS PLGEEMRDRARAHVDALRTHLA PFSDELQRLAARLEALKENGG ARLAEFHAKATEHLSTLSEKAK PALEDLRQGLL PVLESFKVSFLSALEEFTKKLNTQ
Stability & Storage	Lyophilized protein: stable at -20°C for several months. Reconstituted protein: stable at 2-8°C for several days.
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP].

	A Material Safety Data Sheet is provided.
Application	The increased intrinsic fluorescence of MSP1E3D1 "Ultra Bright" allows easy detection and monitoring of nanodisc properties without exogenous labels, and enables its use as an intrinsic FRET donor for labeled proteins. It further supports binding studies via analytical ultracentrifugation (AUC) and microscale thermophoresis (MST) without the need for extrinsic fluorophores. Beyond its fluorescence advantages, MSP1E3D1 "Ultra Bright" is suited for Cryo-EM, ELISA, protein-protein interaction studies, and cell-free expression systems, supporting a broad range of biomedical, proteomic, and protein biochemistry applications.
Recommended usage	ELISA: dilute 1:100
Manufacturer	Cube Biotech

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.