

NativeMP™ Wash Buffer (10x)

Cat. No.: 91001

Description	NativeMP™ Wash Buffer (10x) is a 10x concentrated wash buffer suitable for membrane proteins. It can be used in combination with Strep-Tactin®XT, Rho1D4 or anti-DYKDDDDK purification plates or resins.
Components	200 mM HEPES ,1.5 M NaCl, pH 7.5
Amount	100 mL
Form	Solution
Concentration	10x concentrated
Application	Dilute 1 volume of 10x concentrated solution with 9 volumes of demineralized or distilled water.
Stability	12 months after shipping
Storage	2-8 °C (for both concentrated and working solution)
Shipping	Room temperature, blue ice or dry ice
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Manufacturer	Cube Biotech

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.