

PureCube Ni-NTA FPLC column

Cat. No.: 31302, 31306

Description	PureCube Ni-NTA FPLC columns contain PureCube Ni-NTA Agarose. A high density of nickel ions (Ni) loaded onto the ligand nitrilotriacetic acid (NTA) ensures a high binding capacity. The columns are applicable for purification of proteins fused to a His-tag via chromatography workstations using 10-32 connections, syringes, or peristaltic pumps.			
Support	7.5% agarose, highly cross-linked			
Form	Pre-packed in buffer (20 mM sodium acetate with 20% ethanol, pH 6.6)			
Bead size	40 µm, spherical			
Binding capacity	80 mg protein/ml resin. Binding capacity was determined with a 27 kDa 6xHis-tagged GFP protein. Please note that the binding capacity is protein dependent.			
Compatible reagents	1 mM DTT, 1 mM EDTA, 100% methanol, 100% ethanol, 8 M urea, 6 M guanidinium hydrochloride, 30% (v/v) acetonitrile.			
Application	The PureCube FPLC column is compatible with almost all commonly used FPLC systems like ÄKTA™.			
	FPLC column volume	Physical dimension	Flow rate	Pressure stability
	1 ml (31302)	Inner diameter: 0.7 cm Height: 2.5 cm	0.5-2.0 ml/min (max. 6 ml/ml)	Maximum 5 bar (72 psi)
	5 ml (31306)	Inner diameter: 1.6 cm Height: 2.5 cm	0.5-2.0 ml/min (max. 6 ml/ml)	Maximum 3 bar (42 psi)
Stability	2 years after shipping			
Storage	2-8 °C			
Shipping	Room temperature			
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.			
Manufacturer	Cube Biotech			

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.