

NativeMP™ DYKDDDDK Elution Buffer

Cat. No.: 91030

Description	NativeMP™ DYKDDDDK Elution Buffer contains DYKDDDDK peptide for elution. It is especially suitable for membrane proteins and can be used in combination with anti-DYKDDDDK purification plates.
Components	20 mM HEPES , 150 mM NaCl, 0.5 mg/ml DYKDDDDK peptide, pH 7.5
Amount	25 mL
Form	Solution
Concentration	Ready-to-use
Stability	2 years after shipping
Storage	-20°C for long-term storage; 2-8 °C for up to two weeks
Shipping	Blue ice or dry ice
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Manufacturer	Cube Biotech

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.