

MSP1E3D1 “Dark”, human

Cat. No.: 26362

Description	MSP1E3D1 “Dark” is a membrane scaffold protein (MSP) derived from apolipoprotein A-I that wraps around a phospholipid bilayer patch to form stable, water-soluble nanodiscs of 12–14 nm in diameter, enabling the stabilization of membrane proteins in a defined lipid environment. Its defining feature is the replacement of all tryptophan residues with phenylalanines, resulting in a significantly reduced intrinsic fluorescence and UV absorbance at 280 nm, making accurate concentration determination of MSP-stabilized membrane proteins easier than ever.
Form	Lyophilized
Amount	2 mg
Available origins	Human
Purity	> 90% (Determined by SDS-PAGE)
Number of amino acids	256
Molecular mass	29.86 kDa
Extinction coefficient (in H₂O) ϵ_{280}	Ext. coefficient 10430 Abs 0.1% (=1 g/l) 0.349
Buffer	20 mM TRIS pH 7.4, 100 mM NaCl, 0.5 mM EDTA
Affinity tag	No tag
Sequence	G STFSKLREQLG PVTQEF F DNLEKETEGLRQEMS KDLEEVKAKVQ PYLDDFQKK F QEEMELYRQKVE PLRAELQEGARQKLHELQEKLS PLGEEMRDRARAHVDALRTHLA PYLDDFQKK F QEEMELYRQKVE PLRAELQEGARQKLHELQEKLS PLGEEMRDRARAHVDALRTHLA PYSDELQRQLAARLEALKENG ARLAEYHAKATEHLSTLSEKAK PALEDLRQGLL PVLESFKVSFLSALEEYTKKLNTQ
Stability & Storage	Lyophilized protein: stable at -20°C for 24 months. Reconstituted protein: stable at 2-8°C for several days.
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

Application	MSP1E3D1 "Dark" enables reliable UV-Vis based concentration determination and reduces background interference in fluorescence-based assays and label-free detection methods, offering a clear advantage over standard MSPs and alternative solubilization agents such as AASTY or SMA. It is further suited for Cryo-EM, ELISA, protein-protein interaction studies, and cell-free expression systems, supporting a broad range of biomedical, proteomic, and protein biochemistry applications.
Recommended usage	
Manufacturer	Cube Biotech

For research use only**Trademark information**

The owners of trademarks marked by “®” or “TM” are identified at <http://www.cube-biotech.com/information/patents>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale Cube Biotech grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.cube-biotech.com/information/patents> or upon inquiry at info@cube-biotech.com or at Cube Biotech GmbH, Creative-Campus-Allee 12, D-40789 Monheim, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.