



ORBISCOPE

Powered by Polyscope Polymers

Certificate of Analysis **SMALP 30010EP**

Product information

Product	SMALP 30010EP
Description	Aqueous solution of an esterified SMA® polymer
Lot no.	PIT 093780-LN-A

Analyses result	Value *	Unit
Dry solids content	22.8	% w/w
pH	8.0	
Absolute viscosity	<500	mPa·s
Color		Yellow

Orbiscope certifies that the above information is correct.

Linda Nijssen
Lab Technician

Stefan Scheidelaar, PhD
Technical Product Manager

* Analytical procedures are available upon request



ORBISCOPE

Powered by Polyscope Polymers

Prins de Lignestraat 23 • 6161 CZ Geleen • The Netherlands

T +31 (0)46 - 750 00 10 • info@orbiscope.com

www.orbiscope.com